



be INSPIRED drive DISCOVERY stay GENUINE

version 1.16.5

(i) ABOUT

? HELP

Product Group

Taq DNA Polymerase

Polymerase/Kit

Taq DNA Polymerase with Standard Taq Buffer

Primer Concentration (nM)

200

Primer Pairs

Type/paste sequences below or load primer pairs from file:

Choose File No file chosen

Interleaved data (one primer per line)

Plfwd; AGCGGATAACAATTTCACACAGGA; Plrev; GTAAAACGACGGCCAGT

pBeta; AGCGGATAACAATTTCAC

P3fwd; AGCGGATAAGGGCAATTTCAC; P3rev; GTAAAACGACGGCCA

5 primers processed. Errors: 1 Warnings: 1



ID	Sequence	Tm	Та	Notes
P1fwd	AGCGGATAACAATTTCACACAGGA	58	49	ОК
P1rev	GTAAAACGACGGCCAGT	54	49	ОК
pBeta	AGCGGATAACAATTTCAC	48		Primer 2 missing Both primers need to be longer than 7 nt
P3fwd	AGCGGATAAGGGCAATTTCAC	57	45	OKTm difference is greater than the recommended limit of 5 °C.
P3rev	GTAAAACGACGGCCA	50	45	OKTm difference is greater than the recommended limit of 5 °C.

Tech Support Feedback NEB Overview Site Map

Trademarks Legal and Disclaimers Privacy Cookie Policy Terms of Use

US